

Rat ALK-4/ACVR1B protein (Recombinant) (C-His) (STJP002979)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29381
Gene Symbol	Acvr1b
UniProt ID	ACV1B_RAT
Immunogen Region	Met1-Glu126
Immunogen Sequence	CDLPQTHSLGNRRALILLAQ MGRISHFSLKDRHDFGFPE EEFDGHQFQKAQAISVLHEM IQQTFNLFSTEDSSAAWEQS LLEKFSTELYQQLNDLEACV IQEVGVEETPLMNEDSILAV RKYFQRITLYLTEKKYSPCA WEVVRAEIMRSLSFSTNLQK RLRRKD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Rat ALK-4/ACVR1B Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 40-50 kDa.