

Human KLHL20 protein (Recombinant) (N-His) (STJP002980)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Supplied as a 0.22 µg/ml filtered solution in PBS pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27252
Gene Symbol	KLHL20
UniProt ID	KLH20_HUMAN
Immunogen Region	Met303-Glu605
Immunogen Sequence	VSLRIQEGCSHLPSCCPSKE QEPPEEWLKWSSASVSPPEP LSHTHHAESCRASKDGPLNS RAISPWSYELDRDLNRVPQD LYHARCLCPHCVSLQTGSHM DPLGNSVPLYHNQTVFYRRP CHGEEGTHRRYCLERRLYRV SLACVCVRPRVMA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human KLHL20 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 35-45 kDa.