

Mouse IL-27A protein (Recombinant) (No-Tag) (STJP002984)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	246779
Gene Symbol	IL27
UniProt ID	IL27A_MOUSE
Immunogen Region	Phe29-Ser234
Immunogen Sequence	HPAGKRPCKMQAFRIWDTNQ KTFYLRNNQLIAGYLQGPNT KLEEKIDMVPIDFRNVFLGI HGGKLCCLSCVKSGDDTKLQL EEVNITDLNKNKEEDKRFTF IRSETGPTTSFESLACPGWF LCTTLEADHPVSLTNTPKPEP CTVTKFYFQEDQ

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Mouse IL-27A Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 25-35 kDa.