

Human EREG protein (Recombinant) (C-His) (STJP002999)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2069
Gene Symbol	EREG
UniProt ID	EREG_HUMAN
Immunogen Region	Val63-Leu108
Immunogen Sequence	QGVFEDCCLAYHYPIGWAVL RRAWTYRIQEVSGSCNLPAA IFYLPKRHRKVCGNPKSREV GRAMKLLDARNKVFAKLHHN TQTFQAGPHAVKKLSSGNSK LSSSKFSNPISSSKRNVSLL ISANSGL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human EREG Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 35-40kDa.