

Human Gamma-glutamyl hydrolase/GGH protein (Recombinant) (C-His) (STJP003005)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of 20mM HEPES, 150mM NaCl, 1mM TCEP, pH 8.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8836
Gene Symbol	GGH
UniProt ID	GGH_HUMAN
Immunogen Region	Arg25-Asp318
Immunogen Sequence	PARTAPARVPTLAVPAISLP DDVRRRLKDLERDSLTEKEC VKEKLNLLHEFLQTEIKNQL CDLETCLRKEELSEEGLAK VKSLLNKDLSLENGAHAYNR EVNGRLENGNQARSEARRVG MADANSPPKPLSKPRTPRRS KSDGEAKRSRDPASASQVT GIRAEPPSPRITRKSTRQT TITSHFAKGPARKRPQEESE RAKSDESIKEEDKDQDEKRR RVTSRERVARPLPAEEPER

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human Gamma-glutamyl hydrolase/GGH Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.