

Human ITPRIPL1/KIAA1754L/CD3L1 protein (Recombinant) (No-Tag) (STJP003008)

GENERAL INFORMATION

| | |
|---------------|----------|
| Product Type | Proteins |
| Host / Source | E.coli |

PRODUCT PROPERTIES


| | |
|---------------------|--|
| Formulation | Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4. |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|--------------------|---|
| Gene ID | 150771 |
| Gene Symbol | ITPRIPL1 |
| UniProt ID | IPIL1_HUMAN |
| Immunogen Region | His25-103Gly |
| Immunogen Sequence | PETLCGAEIYDALQFVCGDR GFYFNKPTGYGSSRRAPQT GIVDECCFRSCDLRRLEMYC APLKPAKSA |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|



Recombinant ITPRIPL1/KIAA1754L/CD3L1 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.