

Human CXCL4/PF-4 protein (Recombinant) (C-His) (STJP003017)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 M filtered solution of 10 mM Sodium Citrate, 150 mM NaCl, pH 5.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5196
Gene Symbol	PF4
UniProt ID	PLF4_HUMAN
Immunogen Region	Glu32-Ser101
Immunogen Sequence	APLAADPTACCFSYTSRQI PQNFIADYFETSSQCSKPSV IFLTKRGRQVCADPSEEWVQ KYVSDLELSA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human PF-4/CXCL4/SCYB4 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 10-15 kDa.