


Monkeypox virus D6L protein (Recombinant) (C-His) (STJP003032)

GENERAL INFORMATION		
	Product Type	Proteins
	Host / Source	Pichia
PRODUCT PROPERTIES		
	Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
	Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Immunogen Region	Val21-Lys126
	Immunogen Sequence	ARATNVGRECCLDYFKGAIP IRKLVSWYKTSVECSRDAIV FLTVQGKLCADPKDKHVKK AIRLVKNPRP
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant Monkeypox virus D6L Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 14-15 kDa.