


Mouse IL-20 protein (Recombinant) (C-His) (STJP003035)

GENERAL INFORMATION	
Product Type	Proteins
Host / Source	HEK293 cells
PRODUCT PROPERTIES	
Formulation	Lyophilised from a 0.22 Mu m filtered solution of 20mM PB, 500mM NaCl, pH 8.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	58181
Gene Symbol	IL20
UniProt ID	IL20_MOUSE
Immunogen Region	Leu25-Leu176
Immunogen Sequence	QDSTSDLIPAPPLSKVPLQQ NFQDNQFQGKQWYVVLGNA ILREDKDPQKMYATIELKE DKSYNVTSLVFRKKKCDYWI RTFVPGCQPGFTLGNIKSYPGLTSYLVRVSTNYNQHAM VFFKKVSNREYFKITLYGR TKELTSELKENFIRFSKSLG LPENHIVFPVIDQCIDG
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant Mouse IL-20 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 15-25 kDa.