

Human CD81 protein (Recombinant) (C-His) (STJP003040)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	975
Gene Symbol	CD81
UniProt ID	CD81_HUMAN
Immunogen Region	Phe113-Lys201
Immunogen Sequence	CNLSQTHSLNNRRTLMLMAQ MRRISPFSCCLKDRHDFEFPQ EEFDGNQFQKAQAISVLHEM MQQTFNLFSTKNSSAAWDET LLEKFYIELFQQMNDLEACV IQEVGVEETPLMNEDSILAV KKYFQRITLYLMEKKYSPCA WEVVRAEIMRSLSFSTNLQK RLRRKD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human CD81 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 35-45 kDa.