

Mouse IL-19 protein (Recombinant) (C-His) (STJP003041)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	329244
Gene Symbol	IL19
UniProt ID	IL19_MOUSE
Immunogen Region	Leu25-Ala176
Immunogen Sequence	CDLPQTHSLGNRRALILLAQ MGRISPFSLKDRHDFGFPO EEFDGNQFQKAQAISVLHEM IQQTFNLFSTKDSSATWEQS LLEKFSTELNQLNDLEACV IQEVGVEETPLMNVDSILAV KKYFQRITLYLTEKKYSPCA WEVVRAEIMRSFSLSKIFQE RLRRKE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Mouse IL-19 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 25-40 kDa.