

Mouse CCL17/TARC protein (Recombinant) (C-His) (STJP003042)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.22 M µ filtered solution in PBS (pH 7.4). Normally 8% Trehalose is added as protectant before lyophilisation.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	20295
Gene Symbol	Ccl17
UniProt ID	CCL17_MOUSE
Immunogen Region	Ala24-Pro93
Immunogen Sequence	LGCDLPQTHSLSNRRTLMM AQMGRISPFSCCLKDRHDFGF PQEEFDGNQFQKAQAISVLH EMIQQTFFNLSTKDSSATWD ETLLDKFYTELYQQNDLEA CMMQEVGVEDTPLMNVDSIL TVRKYFQRITLYLTEKKYSP CAWEVVRAEIMRSFSLSANL QERLRRKE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Mouse CCL17/TARC Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 10-15 kDa.