

Human AMHR2/MISR2 protein (Recombinant) (C-His) (STJP003057)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	269
Gene Symbol	AMHR2
UniProt ID	AMHR2_HUMAN
Immunogen Region	Pro18-Ser144
Immunogen Sequence	MQLKPMEINPEMLNKNVLSRL GVAGQWRFVDVLGLEEESLG SVPAPACALLLLFPLTAQHE NFRKKQIEELKGQEVSPKVY FMKQTIGNSCGTIGLIHAVA NNQDKLGFEDGSVLKQFLSE TEKMSPEDRAKCFEKNEAIQ AAHDAVAQEGQCRVDDKVN FHFILFNNVDGHLYELDGRMP FPNVHGASSEDTLKDAKV CREFTEREQGEVRFSAVALC KAA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human AMHR2/MISR2 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.