

Mouse Leptin/LEP protein (Recombinant) (C-His) (STJP003060)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 M filtered solution of 20mM PB, 300mM NaCl, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	16846
Gene Symbol	Lep
UniProt ID	LEP_MOUSE
Immunogen Region	Val22-Cys167
Immunogen Sequence	SLGCDLPQTHSLGNRRALIL LAQMRRISPFSLKDRHDFE FPQEEFDDKQFQKAQISVL HEMIQQTFLNLFSTKDSSAAL DETLLDEFYIELDQQLNDLE SCVMQEVGVIESPLMYEDSI LAVRKYFQRITLYLTEKKYS SCAWEVVRAEIMRSFSL SIN LQKRLKSKE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Mouse Leptin/LEP Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 15-25 kDa.