


Human Macrophage migration inhibitory factor/MIF protein (Recombinant) (No tag) (STJP003061)

GENERAL INFORMATION		
	Product Type	Proteins
	Host / Source	E.coli
PRODUCT PROPERTIES		
	Formulation	Lyophilised from a 0.22 Mu m filtered solution of 20mM PB, 150mM NaCl, pH 7.4/20mM PB, 150mM NaCl, pH 7.4, 50% Glycerol
	Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	4282
	Gene Symbol	MIF
	UniProt ID	MIF_HUMAN
	Immunogen Region	Met1-Ala115
	Immunogen Sequence	PGPPAGSPRVSSDPRADLDS AVLLTRSLADTRQLAAQMR DKFPADGDHSLDSLPTLAMS AGTLGSLQLPGVLTRLRVDL MSYLRHVQWLRRAGGPSLKT LEPELGALQARLERLLRRLQ LLMSRLALPQAAPDQPVIP L GPPASAWGSIRAAHAILGGL HLTLDWAVRGLLLLKTRL
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant Human Macrophage migration inhibitory factor/MIF Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.