

Mouse Midkine/MDK protein (Recombinant) (C-His) (STJP003066)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	17242
Gene Symbol	Mdk
UniProt ID	MK_MOUSE
Immunogen Region	Lys23-Asp140
Immunogen Sequence	MKLLVLAVLLTVAAADSGIS PRAVWQFRKMIKCVIPGSDP FLEYNNGCYCGLGGSGTPV DELDKCCQTHDNCYDQAKKL DSCKFLLDNPYTHYYSYSCS GSAITCSSKNKECEAFICNC DRNAAICFSKAPYNKAHKNL DTKKYCQS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Mouse Midkine Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 10-15 kD.