

## Human Ubiquitin thioesterase L1/UCHL1 protein (Recombinant) (No tag) (STJP003067)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Host / Source | E.coli   |

### PRODUCT PROPERTIES


|                     |                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------|
| Formulation         | Lyophilised from a 0.22 Mu m filtered solution of 20mM Tris-HCl, 250mM NaCl, 0.05% Tween80, 1mM DTT, pH 8.5. |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.               |

### TARGET INFORMATION

|                    |                                                                                                                                                                                                                                           |
|--------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">7345</a>                                                                                                                                                                                                                      |
| Gene Symbol        | <a href="#">UCHL1</a>                                                                                                                                                                                                                     |
| UniProt ID         | <a href="#">UCHL1_HUMAN</a>                                                                                                                                                                                                               |
| Immunogen Region   | Met1-Ala223                                                                                                                                                                                                                               |
| Immunogen Sequence | MSPILGYWKIKGLVQPTRLL LEYLEEKYEEHLYERDEGDK WRNKKFELGLEFPNLPYYID<br>GDVKLTQSMAIRYIADKHN MLGGCPKERAISMLEGAVL DIRYGVSRAYSKDFETLKV<br>DFLSKLPEMLKMFEDRLCHK TYLNGDHVTHPDFMLYDALD VVLYMDPMCLDAFPKLVCFK<br>KRIEAIQIDKYLKSSKYIA WPLQGWQATFGGGDHPPK |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|



Recombinant Human Ubiquitin thioesterase L1/UCHL1 protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 25-35 kDa.