

Human Follistatin-like 1/FSTL1 protein (Recombinant) (No-Tag) (STJP003068)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% Trehalose is added as protectant before lyophilisation.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	11167
Gene Symbol	FSTL1
UniProt ID	FSTL1_HUMAN
Immunogen Region	Glu21-Ile308
Immunogen Sequence	MSDSEVNQEAKPEVKPEVKP ETHINLKVSDGSSEIFFKIK KTTPLRRLMEAFKRQGGKEM DSLRFLYDGI RIQADQTPED LDMEDNDIIEAHREQIGG

ADDITIONAL INFORMATION


Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human Follistatin-like protein 1/FSTL1 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 45-60kDa.



Recombinant Human Follistatin-like protein 1/FSTL1 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 45-60kDa.



The purity of Human FSTL1 is greater than 95% as determined by SEC-HPLC.