

Mouse FGF-1/aFGF protein (Recombinant) (N-His (IL-23a) &No tag (IL-12b)) (STJP003074)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.2 M filtered solution of 20 mM Tris-HCl, pH 7.5, 50 mM Na ₂ SO ₄ , 0.5 mM DTT, 1 mM EDTA.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	14164
Gene Symbol	Fgf1
UniProt ID	FGF1_MOUSE
Immunogen Region	Phe16-Asp155
Immunogen Sequence	MWELEKDVYVVEVDWTPDAP GETVNLTCDTPEEDDITWTS DQRHGVIGSGKTLTITVKEF LDAGQYTCHKGGETLSHSHL LLHKKENGIWSTEILKNFKN KTFLEAPNYSGRFTCSWL VQRNMDLKFNIKSSSSSPDS RAVTCGMASLSAEKVTLDQR DYKYSVSCQEDVTCPTAEE TLPIELALEARQQNKYENYS TSFFIRDIKPDPPKNLQMK PLKNSQVEVSWEYPDSWST

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Mouse FGF-1/aFGF Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 15-25 kDa.