

Rat IL-10 protein (Recombinant) (C-His) (STJP003075)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25325
Gene Symbol	Il10
UniProt ID	IL10_RAT
Immunogen Region	Ser19-Asn178
Immunogen Sequence	WPESPKAQASSVPTAQPAE GSLAKATTAPATTRNTGRGG EEKKKEKEKEEQEERETKTP ECPSHTQPLGVYLLTPAVQD LWLRDKATFTCFVVGSDLKD AHLTWEVAGKVPTGGVEEGL LERHSNGSQSQHSRLTPRS LWNAGTSVTCTLNHPSLPPQ RLMALREPAQAQAPVKLSNL LASSDPPEAASWLLCEVSGF SPPNILLMWLEDQREVNTSG FAPARPPPQPRSTTFWAWS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Rat IL-10 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.