

Schistosoma japonicum GST protein (Recombinant) (N-His) (STJP003077)

GENERAL INFORMATION

| | |
|---------------|--------------|
| Product Type | Proteins |
| Host / Source | HEK293 cells |

PRODUCT PROPERTIES


| | |
|---------------------|------------------------------------------------------------------------------------------------|
| Formulation | Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4. |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunogen Region | Met1-Lys218 |
| Immunogen Sequence | DTAIRVFAIPPSFASIFLTK STKLTCLVTDLTTYDSVTIS WTRQNGEAVKTHTNISESHP NATFSAVGEASICEDDWNSG ERFCTVTHDLPSPKQTI SRPKGVALHRPDVYLLPPAR EQLNLRRESATITCLVTGFSP ADVFVQWMQRGQPLSPEKYV TSAPMPEQPAPGRYFAHSIL TVSEEEWNTGETYTCVVAHE ALPNRVTERTVDKSTGKPTL YNVSLVMSDTAGTCY |

ADDITIONAL INFORMATION

| | |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|



Recombinant Schistosoma japonicum GST Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 25-35 kDa.