

Human IgA1 protein (Recombinant) (C-His) (STJP003080)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% Trehalose is added as protectant before lyophilisation.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	Cys122-Tyr353
Immunogen Sequence	LDKYGVSDYHKNLINNAKTVEGVKDLQAQVVESAKKARISEATDGLSDFLKSQTPAEDTV KSIELAEAKVLNRELDKYGVSDYYKNLINNAKTVEGVKALIDEILAAPKTDITYKLILN GKTLKGETTTEAVDAATAEKVFKQYANDNGVDGEWYDDATKTFTVTEKPEVIDASELTP AVTTYKLVLINGKTLK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant Human IgA1 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 30-40 kDa.