

Human Lactate Dehydrogenase A/LDHA protein (Recombinant) (C-His) (STJP003082)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Supplied as a 0.22 µg/ml filtered solution in 20mM Tris, 250mM NaCl, 10% Glycerol, pH 8.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3939
Gene Symbol	LDHA
UniProt ID	LDHA_HUMAN
Immunogen Region	Met1-Phe332
Immunogen Sequence	QSYSIYENYGNKEFNATHA AFEYSKISIGGTPALDRRVQD VNDTISDVKQKWRCVVYPGN GFVSASIFGFQAEVGPNNTR SIRKFNTMRQCIDFTFSDVI NIDIYNPCAPNINNTecQF LKSVL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human Lactate Dehydrogenase A/LDHA Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 35-45 kDa.