

Human IgD protein (Recombinant) (C-His) (STJP003085)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% Trehalose is added as protectant before lyophilisation.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	Arg99-Lys384
Immunogen Sequence	QQTRQEADRG CETLVVQHG H CSYTFLLPKSEPCPPGPEVS RDSNTLQRESLANPLHLGKL PTQQVKQLEQALQNNTQWLK KLERAIKTLRSKLEQVQQQ MAQNQTAPMLELGTSLNQT TAQIRKLTDMEAQLLNQTSR MDAQMPETFLSTNKLENQLL LQRQKLQQLQGQNSALEKRL QALETKQQEELASILSKKAK LLNTLRSQSAALTNIERGLR GVRHNSLLQDQQHSLRQL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human IgD Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 45-60kDa.