

Human TNFSF13/APRIL/CD256 protein (Recombinant) (No-Tag) (STJP003096)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8741
Gene Symbol	TNFSF13
UniProt ID	TNF13_HUMAN
Immunogen Region	Lys111-leu250
Immunogen Sequence	SAGATAADGPCALRELSVDL RAERSVLIPETYQANNCQGV CGWPQSDRNPRYGNHVLLL KMQVRGAALARPPCCVPTAY AGKLLISLSEERISAHVVPN MVATECGCR

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human TNFSF13/APRIL/CD256 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 50 kDa.