

Human CXC-chemokine receptor 1/CXCR1 protein (Recombinant) (C-hFc) (STJP003103)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	Human

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% Trehalose is added as protectant before lyophilisation.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3577
Gene Symbol	CXCR1
UniProt ID	CXCR1_HUMAN
Immunogen Region	Met1-Lys39
Immunogen Sequence	MSNITDPQMWFDDLNFTGM PPADEDYSPCMLETETLNK

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human CXC-chemokine receptor 1/CXCR1 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 45-50 kDa.