

Human Mature MIS/AMH protein (Recombinant) (No-Tag) (STJP003106)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	Human

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of 50mM acetic acid, pH 3.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	268
Gene Symbol	AMH
UniProt ID	MIS_HUMAN
Immunogen Region	Leu25-Arg560
Immunogen Sequence	SAGATAADGPCALRELSVDL RAERSVLIPETYQANNCQGV CGWPQSDRNPRYGNHVLLL KMQVRGAALARPPCCVPTAY AGKLLISLSEERISAHHPVN MVATECGCR

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---