

Human CTHRC1 protein (Recombinant) (C-His) (STJP003107)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	Human

PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	115908
Gene Symbol	CTHRC1
UniProt ID	CTHR1_HUMAN
Immunogen Region	Ser31-Lys243
Immunogen Sequence	SEIPKGKQKAQLRQREVVDL YNGMCLQGPAGVPGRDGSPG ANGIPGTPGIPGRDGFKEK GECLRESFEESWTPNYKQCS WSSLNYGIDLGKIAECTFTK MRSNSALRVLFSGSLRLKCR NACCQRWYFTFNGAECGSL PIEAIIYLDQGSPEMNSTIN IHRTSSVEGLCEGIGAGLVD VAIWVGTCSDYPKGDASTGW NSVSRIIEELPK

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human CTHRC1 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 25-35kDa.