

Human Lymphotactin/XCL1 protein (Recombinant) (C-6*His) (STJP003108)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	Human

PRODUCT PROPERTIES


Formulation	Lyophilised from 0.22 Mu m filtered solution in PBS (pH 7.4). Normally 8% Trehalose is added as protectant before lyophilisation.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6375
Gene Symbol	XCL1
UniProt ID	XCL1_HUMAN
Immunogen Region	Val22-Gly114
Immunogen Sequence	VGSEVSDKRTC VSLTTQRLP VSRIKTYTITEGSLRAVIFI TKRGLKVCADPQATWVRDVV RSMDRKSNTNRNNMIQTKPTG TQQSTNTAVTLTG

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human Lymphotactin/XCL1 Protein was determined by SDS-PAGE with Coomassie Blue, showing a band at 20-25kDa.