

## Human Caspase-7/CASP7 protein (Recombinant) (C-His) (STJP003111)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	Human

### PRODUCT PROPERTIES


Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">840</a>
Gene Symbol	<a href="#">CASP7</a>
UniProt ID	<a href="#">CASP7_HUMAN</a>
Immunogen Region	Met1-Gln303
Immunogen Sequence	MADDQGCIEEQGVDSANED SVDAPDRSSFVPSLFSSKKK KNVMTMRSIKTTRDRVPTYQY NMNFEKLGKCIINNKNFDFK VTGMGVRNGTDKDAEALFKC FRSLGFDVIVYNDSCAKMQ DLLKKASEEDHTNAACFACI LLSHGEEENVYKDGVTPIK DLTAHFRGDRCKTLLEKPKL FFIQACRGTELDGDIQADSG PINDTDANPRYKIPVEADFL FAYSTVPGYYSWRSPGRGS

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant Human Caspase-7/CASP7 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.