

Human CD232/PLXNC1 protein (Recombinant) (N-His) (STJP003313)

GENERAL INFORMATION

| | |
|---------------|-----------------------------|
| Product Type | Proteins |
| Applications | ELISA/Immunogen/SDS-PAGE/WB |
| Host / Source | E.coli |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Formulation | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol. |
| Molecular Weight | 38.51 kDa |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |

TARGET INFORMATION

| | |
|------------------|-----------------------------|
| Gene ID | 10154 |
| Gene Symbol | PLXNC1 |
| UniProt ID | PLXC1_HUMAN |
| Immunogen Region | Ala35-Thr370 |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|