

Human CGB3 protein (Recombinant) (N-His) (STJP003936)

GENERAL INFORMATION

| | |
|---------------|-----------------------------|
| Product Type | Proteins |
| Applications | ELISA/Immunogen/SDS-PAGE/WB |
| Host / Source | E.coli |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Formulation | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol. |
| Molecular Weight | 16.62 kDa |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |

TARGET INFORMATION

| | |
|------------------|----------------------------------|
| Gene ID | 1082/93659/94115 |
| Gene Symbol | CGB3.CGB5.CGB8 |
| UniProt ID | CGB3_HUMAN |
| Immunogen Region | Pro31-Gln165 |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|