

Human SLC16A1 protein (Recombinant) (N-His) (STJP004921)

GENERAL INFORMATION

| | |
|---------------|-----------------------------|
| Product Type | Proteins |
| Applications | ELISA/Immunogen/SDS-PAGE/WB |
| Host / Source | E.coli |

PRODUCT PROPERTIES

| | |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formulation | Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol. |
| Molecular Weight | 8.61 kDa |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |

TARGET INFORMATION

| | |
|------------------|----------------------------|
| Gene ID | 6566 |
| Gene Symbol | SLC16A1 |
| UniProt ID | MOT1_HUMAN |
| Immunogen Region | Asn444-Val500 |

ADDITIONAL INFORMATION

| | |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|