

Human CRMP2/DPYSL2 protein (Recombinant) (N-His) (STJP005490)

GENERAL INFORMATION	
Product Type	Proteins
Applications	ELISA/Immunogen/SDS-PAGE/WB
Host / Source	E.coli
PRODUCT PROPERTIES	
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Molecular Weight	54.66 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
TARGET INFORMATION	
Gene ID	1808
Gene Symbol	DPYSL2
UniProt ID	DPYL2_HUMAN
Immunogen Region	Thr13-Glu490
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.