

## Human IL34 protein (Recombinant) (C-His) (STJP005611)

### GENERAL INFORMATION

|               |                             |
|---------------|-----------------------------|
| Product Type  | Proteins                    |
| Applications  | ELISA/Immunogen/SDS-PAGE/WB |
| Host / Source | Mammalian Cells             |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Formulation         | Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.   |
| Molecular Weight    | 28.25 kDa   |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |

### TARGET INFORMATION

|                  |                            |
|------------------|----------------------------|
| Gene ID          | <a href="#">146433</a>     |
| Gene Symbol      | <a href="#">IL34</a>       |
| UniProt ID       | <a href="#">IL34_HUMAN</a> |
| Immunogen Region | Asn21-Pro242               |

### ADDITIONAL INFORMATION

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|