

Human CD85d/LILRB2 protein (Recombinant) (N-His) (STJP006684)

GENERAL INFORMATION	
Product Type	Proteins
Applications	ELISA/Immunogen/SDS-PAGE/WB
Host / Source	E.coli
PRODUCT PROPERTIES	
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Molecular Weight	24.92 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
TARGET INFORMATION	
Gene ID	10288
Gene Symbol	LILRB2
UniProt ID	LIRB2_HUMAN
Immunogen Region	Leu53-Val255
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.