

3C Protease protein (Recombinant) (STJP006744)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	50mM Tris-HCl, pH 8.0, 150mM NaCl, 10mM EDTA, 1mM DTT, 50% glycerol
Molecular Weight	47.46 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week. Store at -20 to -80 °C for twelve months from the date of receipt.

TARGET INFORMATION

Immunogen Region	3C Protease is cloned from Human rhinovirus and expressed in E.coli.
------------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--