

Human MDA5/IFIH1 protein (Recombinant) (N-His) (STJP006748)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/Immunogen/SDS-PAGE/WB
Host / Source	Insect Cells

PRODUCT PROPERTIES

Formulation	Lyophilized from a solution in 20mM Tris-HCl, pH 8.0, 100mM NaCl, 1mM DTT, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Molecular Weight	118.49 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID	64135
Gene Symbol	IFIH1
UniProt ID	IFIH1_HUMAN
Immunogen Region	Met1-Asp1025

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--