

Human CADM-140/MDA5 protein (Recombinant) (STJP006772)

GENERAL INFORMATION

Product Type Proteins

PRODUCT PROPERTIES

Formulation 20mM Tris-HCl, pH 8.0, 100mM NaCl, 1mM DTT, 10% Glycerol.
Storage Instruction Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID [64135](#)
Gene Symbol [IFIH1](#)
UniProt ID [IFIH1_HUMAN](#)

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.