

Human DSG1 protein (Recombinant) (STJP006773)

GENERAL INFORMATION

Product Type Proteins

PRODUCT PROPERTIES

Formulation Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.

Storage Instruction Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID [1828](#)

Gene Symbol [DSG1](#)

UniProt ID [DSG1_HUMAN](#)

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.