

Human Calreticulin/CALR protein (Recombinant) (STJP006776)

GENERAL INFORMATION

Product Type Proteins

PRODUCT PROPERTIES

Formulation 0.01M PBS, pH 7.4.

Storage Instruction Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID [811](#)

Gene Symbol [CALR](#)

UniProt ID [CALR_HUMAN](#)

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.