

Human BP180/NC16a protein (Recombinant) (STJP006778)

GENERAL INFORMATION

Product Type Proteins

PRODUCT PROPERTIES

Formulation Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.

Storage Instruction Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Gene ID [1308](#)

Gene Symbol [COL17A1](#)

UniProt ID [COHA1_HUMAN](#)

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.