

CSFV N-Pro protein (Recombinant) (C-His) (STJP006833)

GENERAL INFORMATION		
	Product Type	Proteins
	Applications	ELISA/Immunogen/SDS-PAGE/WB
	Host / Source	Insect Cells
PRODUCT PROPERTIES		
	Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
	Molecular Weight	20.01kDa
	Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
TARGET INFORMATION		
	Immunogen Region	Met1-Ser167
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.