

## HBV C/Capsid protein protein (Recombinant) (C-His) (STJP006849)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/WB/ELISA/Immunogen
Host / Source	E.coli

### PRODUCT PROPERTIES

Formulation	0.01M PBS, pH 7.4.
Molecular Weight	17.92 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 2 to 8 °C for one week.Store at-20 to-80 °C for twelve months from the date of receipt.

### TARGET INFORMATION

Immunogen Region	Met1-Val149
------------------	-------------

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--