

Beta-Hexosaminidase B (HEXB) protein (Recombinant) (STJP006852)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	Mammalian

PRODUCT PROPERTIES

Formulation	PBS pH 7.4, 50% glycerin
Molecular Weight	62 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 2 to 8 °C for one week.Store at-20 to-80 °C for twelve months from the date of receipt.

TARGET INFORMATION

Immunogen Region	Beta-Hexosaminidase B (HEXB) is expressed in Mammalian.
------------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--