

Benzonase Nuclease protein (Recombinant) (STJP006947)

GENERAL INFORMATION

| | |
|---------------|----------|
| Product Type | Proteins |
| Host / Source | E.coli |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Formulation | 25 mM Tris pH 8.0, 100mM NaCl, 50% Glycerol. |
| Molecular Weight | 30.16 kDa |
| Storage Instruction | Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 2 to 8 °C for one week.Store at-20 to-80 °C for twelve months from the date of receipt. |

TARGET INFORMATION

| | |
|------------------|---|
| Immunogen Region | Benzonase Nuclease is cloned from <i>Serratia marcescens</i> and expressed in E.coli. |
|------------------|---|

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|