

Endo S2/Endo-beta-N-acetylglucosaminidase protein (Recombinant) (STJP006959)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from PBS pH7.4, 1mM EDTA, 4%trehalose, 1% mannitol.
Molecular Weight	114.67 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles.Store at 2 to 8 °C for one week.Store at-20 to-80 °C for twelve months from the date of receipt.

TARGET INFORMATION

Immunogen Region	Endo S2/Endo-beta-N-acetylglucosaminidase is cloned from Streptococcus pyogenes and expressed in E.coli.
------------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--