

BVDV Glycoprotein E2 protein (Recombinant) (N-His) (STJP007178)

GENERAL INFORMATION		
	Product Type	Proteins
	Applications	ELISA/Immunogen/SDS-PAGE/WB
	Host / Source	E.coli
PRODUCT PROPERTIES		
	Formulation	0.01M PBS, pH 7.4.
	Molecular Weight	20.34 kDa
	Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
TARGET INFORMATION		
	Immunogen Region	His693-Ser851
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.