

DENV-2 Envelope protein E protein (Recombinant) (N-His-GFP) (STJP007219)

GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/Immunogen/SDS-PAGE/WB
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilized from PBS pH7.4, 0.02%NLS, 1mM EDTA, 4%trehalose, 1% mannitol.
Molecular Weight	77.29 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

TARGET INFORMATION

Immunogen Region	Met281-Ala722
------------------	---------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--