

## Human FGF19 protein (Recombinant) (N-His tag) (STJP007235)

### GENERAL INFORMATION

Product Type	Proteins
Applications	ELISA/Immunogen/SDS-PAGE/WB
Host / Source	E.coli

### PRODUCT PROPERTIES

Formulation	0.01M PBS, pH 7.4.
Molecular Weight	23.54 kDa
Storage Instruction	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

### TARGET INFORMATION

Gene ID	<a href="#">9965</a>
Gene Symbol	<a href="#">FGF19</a>
UniProt ID	<a href="#">FGF19_HUMAN</a>
Immunogen Region	Phe27-Lys216

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--